## Determination of Total Volatile Organic Compounds in Recycled Plastic by Headspace/Gas Chromatography-Mass Spectrometry

ShanhaoXu<sup>1</sup>, ChuanLuo<sup>1</sup>, HaoWang<sup>2</sup>, YinZhang<sup>3</sup>, ChaoZhang<sup>3</sup> and LifengYuan<sup>2\*</sup>

<sup>1</sup> Ningbo Customs District, QinYi Road 66#, NingBo, 315012, P. R. China

<sup>2</sup> Ningbo academy of science and technology for inspection and quarantine, QinYi Road 66#, NingBo, 315012, P. R. China

<sup>3</sup> Joysun Testing service, QinYi Road 66#, NingBo, 315012, P.R. China

**Abstract.** A method of headspace combination with gas chromatography-mass spectrometry was developed for the determination of total volatile organic compounds (TVOC) from recycled plastic. The factors such as equilibrium temperature, equilibrium time and sample uniformity were investigated. Toluene was selected as reference material. The results showed that the optimum conditions were as follows: equilibrium temperature was 80°Cand equilibrium time was 25min. The calibration curves showed good linearity and the correlative coefficients was 0.9975. The method detection limit(S/N=3) was 0.01mg/kg, and the recoveries ranged from 83.5% to 108.3%, with RSDs <6%. The method has the advantages of no pretreatment, accurate results and high sensitivity etc., and can be used for the determination of TVOC from recycled plastic.

### **1** Introduction

Volatile organic compounds (VOCs) refer to organic compounds that are prone to volatilize at room temperature, including alkanes, alkenes, aromatics, halogenated hydrocarbons, alcohols, aldehydes, ketones and so on. TVOC is a term used to represent the pool of volatile organic compounds. Some volatile halogenated compounds are toxic to human body even at low concentration <sup>[1]</sup>, and some VOCs that are hard to degrade can cause bioaccumulation effect, affecting biological system and environment on a large scale <sup>[2]</sup>. Therefore, there are strict control standards and policies for TVOC emission at home and abroad [3]. At present, the detection of VOCs is performed mainly via gas [4] analysis/ chromatography thermal gas chromatography-mass spectrometry <sup>[5]</sup>, headspace-gas chromatography<sup>[6]</sup>.

Recycled plastics generally contain VOCs because of their complex sources, incomplete cleaning, as well as cracking and decomposition during processing. With the ban on imported solid waste and tightening environmental protection policies, higher demands have been proposed for the quality of recycled plastics <sup>[7,8]</sup>. At present, there have been few reports on the detecting methods of VOCs in recycled plastics. The headspace sample (HS) combined with gas chromatography-mass spectrometry (GC-MS) is developed in recent years to detect volatile substances, which was free from tedious pretreatment, interference of solvent impurities and less pollution. This technology has the characteristics of avoiding tedious pretreatment, avoiding interference of solvent impurities, and reducing pollution to the injection port and chromatographic column.

This paper tested TVOC in recycled plastics via headspace-gas chromatography-mass spectrometry that skipped pretreatment and featured high accuracy and sensitivity.

## 2 Experimental part

#### 2.1. Instruments and conditions

Instrument: Gas chromatography-mass spectrometry, GC2030-QP2020NX from SHIMADZU of Japan, with EI source; Headspace sampler: AOC-6000 from SHIMADZU of Japan; Analytical balance, METTLER TOLEDO from TOLEDO of Switzerland, with a sensitivity of 0.1 mg;22 mL headspace bottle; 10*u*L microinjector.

#### 2.2. Reagents and materials

Helium, purity  $\geq$  99.999%; Methanol (HPLC grade from TEDIA of USA); Toluene standard solution (99.7%, Dr.Ehrenstorfer of Germany); Preparation of standard stock solution of toluene (1 mg/mL): 0.1g of toluene standard was weighed and adjusted to 100 mL with methanol;

Preparation of standard solution of toluene  $(20\mu g/mL)$ : 1.0 mL standard stock solution of toluene was removed and adjusted to 50 mL with methanol.

<sup>\*</sup>Corresponding author's e-mail: 18925141@qq. com

<sup>©</sup> The Authors, published by EDP Sciences. This is an open access article distributed under the terms of the Creative Commons Attribution License 4.0 (http://creativecommons.org/licenses/by/4.0/).

analysis via headspace gas chromatography mass

TVOC refers to any organic compound whose boiling

point or initial boiling point is lower than or equal to

250 °C at the standard atmospheric pressure. The

National Standard GB 18581-2009 Limits of Harmful

Substances in Solvent-based Woodware Coatings for

Interior Decoration Materials defines diethyl adipate

(boiling point of 251 °C) as a marker. That is, classify

substances released before the marker as VOC and those after the marker as non-VOC. The sum of substances

released before the marker adds up to the total volatile

the diethyl adipate (boiling point of 251°C) as a marker, and the sum of substances before the marker as TVOC.

This paper followed the above standard and defined

3 Results and discussions

3.1. Selection of markers

organic compounds (TVOC).

spectrometry.

#### 2.3. Methods and conditions

## 2.3.1 Gas chromatography-mass spectrometry conditions

Chromatographic column: DB-624,60 m × 0.25 mm × 1.40  $\mu$ m; Temperature: inlet at 250 °C, auxiliary heating zone at 250 °C, ion source at 200 °C. Keep the state for 4min before rising to 200 °C at the rate of 6 °C/min, and stay at the temperature for 5 min. Next, rise to 240 °C at the rate of 20 °C/min and stay there for 1min. Split sampling in the ratio of 10:1; Carrier gas: Helium at the flow rate of 1.0 mL/min. Ionization: EI (70 ev); Quality scanning: 35 amu-350 amu; Solvent delay: 4min; Scanning mode: Full scan (SCAN).

#### 2.3.2 Headspace conditions:

Balance point temperature: 80 °C; Duration: 25 min; Temperature of sampling needle: 80 °C.

#### 2.3.3 Procedure

0.5 g of sample (to an accuracy of 0.1 mg) was weighed, placed in the headspace bottles and sealded instantly for



Fig.1 Total ion chromatogram of Diethyl Adipate

it is almost impossible to calculate TVOC by accurately quantifying VOCs in all recycled plastics.

#### 3.2. Calculation of TVOC

There are various kinds of compounds in TVOC, which vary with materials, sources and procedures. Therefore,

Table1 VOCs of different kinds of samples

No.	Sample type	TVOC
1	PE	Acetic acid, Decane, Limonene
2	PP	Hexadecane, Nonylene, Toluene, Nonane
3	PA	Aniline, Toluene, Styrene
4	PC	Dichloromethane, Phenol, Styrene, Xylene, Toluene
5	ABS	Styrene, Acrylonitrile, Toluene
6	PS	Dodecane, Trichloroethylene, Styrene, Benzophenone
7	PVC	Hexadecane, Isooctanol, Heptanone
8	PET	Dodecane, Benzothiazole, Toluene

Based on GB 18581-2009 Limits of Harmful Substances in Solvent-based Wood Coatings for Interior Decoration Materials, this paper set the response factor of all VOC substances relative to toluene at 1.0. That is, all VOC substances were regarded as toluene.

#### 3.3. Selection of equilibrium temperature

Volatile organic compound refers to organic compound that is volatile at room temperature, whose concentration in air is closely related to temperature. Generally, the higher the temperature, the higher the concentration. Given the fact that the temperature of recycled plastics peaks at about 60 °C in routine transportation, storage and use, this paper set the equilibrium temperature at 80°C after taking into account analyzing methods of volatile organic compounds in other plastic products <sup>[9,10]</sup>.

#### 3.4. Determining balancing duration

This paper investigated the effects of balancing duration of 5, 10, 15, 20, 25 and 30 min on the analysis results, which showed that the concentration of VOCs in gas phase increased with the extension of equilibrium time. As the sample was in equilibrium at 25min, this paper set the parameter at 25min for time efficiency.



Fig.2 Effects of different equilibrium times

#### 3.5. Investigation of sample homogeneity

Generally, recycled plastic particles are not homogeneous because of their complex sources and components. Therefore, it is necessary to investigate their homogeneity. This paper selected a sample of about 2 kg. After mixing evenly, 8 segments were selected from different parts of the sample to test the proportion of volatile organic compounds via full scanning. This paper also investigated the relative standard deviations of 6 highly-concentrated volatile organic compounds. The results (Table 2) showed that the relative standard deviations of 6 substances were all less than 5%, featuring good uniformity.

Table2 The results of the homogeneity testing (n=8)

No.	compound	retention time /min	RSD(%)
1	Hexane	9.365	1.7
2	Benzene	12.490	2.1
3	Toluene	16.980	3.0
4	Ethylbenzene	20.865	4.9
5	P-xylene	21.155	4.1
6	Isopropyl benzene	23.200	4.6

#### 3.6. Linear relationship and detection limit

Add 1.0, 2.0, 5.0, 10, 20, 50, and 100  $\mu$ L toluene standard solution into seven headspace bottles, seal the bottles immediately and analyze them in descending order of concentration. The peak surface of total ion flow

relative to the concentration of toluene was used to obtain the curve of total volatile organic compounds (Figure 3), with a linear correlation coefficient of 0.9975. Also, the detection limit was determined at 0.01 mg/ kg (MLOD, S/ N = 3).



Fig.3 Linear equation of TVOC

# 3.7. The recovery rate and precision of the method

One PP and one PVC regenerated particle sample were selected for six parallel experiments and standard recovery experiments, the results of which are shown in Table 3. It showed that the relative standard deviation (RSD) of TVOC was less than 6%, with an average recovery rate of 83.5%-108.3% in the range of standard addition. The method features high recovery and precision.

Table 3 Determination results of recoveries and precisions (n=6)

No.	Sample	RSD(%)	Spiked(mg/kg)	Recovery(%)
1	PP	5.62	0.2, 2	89.2, 108.3
2	PVC	3.99	0.2, 2	83.5, 99.8

#### 3.8. Testing of actual samples

Two samples of ABS, PA, PC, PE, PET, PP, PS and PVC were selected and tested using this method. The results are shown in Table 4.

Table4 The results of the samples

No.	Sample type	TVOC(mg/kg)
1	ABS	0.30, 1.21
2	PA	0.56, 0.12
3	PC	2.7, 0.30
4	PE	7.83, 2.13
5	PET	< 0.01, 0.14
6	PP	7.01, 3.91
7	PS	3.75, 0.66
8	PVC	8.82, 3.30

## 4 Conclusion

This paper used headspace-gas chromatography-mass spectrometry to determine TVOC in recycled plastics, with a lower limit of 0.01 mg/ kg, a recovery rate between 83.5% and 108.3%, and a repeatability of less than 6%. The method requires no pretreatment and features accurate results and high sensitivity.

### Acknowledgements

This work was financially supported by the science research project of General Administration of Customs, P. R. CHINA (Fund No.: 2019HK016), the Ningbo public science research project (Fund No.: 2019C50031) and the Zhejiang basis public science research project (Fund No.: LGC20B040001).

#### References

- Spyros K G, Themistokles D L, Anastasia D N. (2001) Comparison of methods for determination of volatile organic compounds in drinking water. Chemosphere, 45: 275-284
- Safarova V I, Sapelnikova S V, Djazhenko E V, . etc. (2004) Gas chromatography-mass spectrometry with headspace for the analysis of volatile organic compounds in waste water. J.Chromatogr. B, 800: 325-330
- 3. KE FY, CHEN Y, HU JY. etc. (2019) Formation mechanism and control technology on volatile organic compounds of regenerated polyester fiber. China Synthetic Fiber Industry, 42(1):49-54
- 4. International Organization for Standardization. (2013) Paints and varnishes Determination of

volatile organic compound (VOC) content — Part 2: Gas-chromatographic method: ISO 11890-2-2013.Published in Switzerland.

- International Organization for Standardization. (2011) Indoor air - Part 6: Determination of volatile organic compounds in indoor and test chamber air by active sampling on Tenax TA sorbent, thermal desorption and gas chromatography using MS or MS-FID: ISO 16000-6-2011. Published in Switzerland.
- International Organization for Standardization. (2005) Paints and varnishes — Determination of the volatile organic compound content of low-VOC emulsion paints (in-can VOC): ISO 17895:2005. Published in Switzerland.
- Luo C., Wang H., Wang Q., etc. (2020) Application of Instrumental in the Identification of Solid Waste Characteristic of Imported recycled ABS Plastic. Earth. Environ.Sci., 508(2020)012205
- Luo C, Dai SJ, Zhang Y., etc. (2021) Identification Method of Solid Waste Characteristic for Imported Recycled Engineering Plastics. Earth Environ. Sci. 621(2021)012050
- GUO XD, CHEN YG, ZHANG PL. etc. (2018) Research on Volatile Organic Compounds in Six Types of Plastic Products. Contemporary Chemical Industry, 47(4):809-813
- JIANG BL, WU XY, XIONG GL. etc. (2019) Research of VOC and Odor of Polyurethane Foams. Polyurethane Industry, 34:20-23